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


AD-8526 Ethernet Converter Instruction Manual

A&D Company, Ltd.

The AD-8526 Ethernet Converter can connect the RS-232C interface of a weighing instrument to the Ethernet (LAN) port of a computer that is not equipped with an RS-232C interface. By using the AD-8526, weight data from a weighing instrument can be managed with a computer connected to the network. The weighing instrument can also be controlled by the connected computer.

Unpacking the AD-8526

AC adapter labels
Affix them to the AC adapter and the connector end of the cable.

IP address label (5 labels)
Write a preset IP address and affix it on the weighing instrument.

IP address protective cover (5 labels)
Affix it on the IP address label.

Note: Only the following cable is included if specified when ordering:
DIN 7-pin to D-Sub 9-pin cable (AX-KO1786-200)
If not specified, the following two communication cables are included:
DIN 7-pin to D-Sub 9-pin cable (AX-KO1786-200), D-Sub 9-pin to D-Sub 9-pin cable (AX-KO2466-200)

Specification

1. Specification		3. RS-232C serial interface	
Operating environment : -10°C to +40°C		Connector: D-Sub 9-pin male	
AC adapter: Confirm that the adapter type is correct for the local voltage and power receptacle type		Transmission form: Asynchronous, bi-directional	
Power consumption: Approx. 11VA (supplied to the AC adapter)		Data format	Baud rate: 600, 1200, 2400*, 4800, 9600, 19200bps
Dimensions: 113(W) x 60(D) x 38(H) mm			Data bits: 7* or 8 bits
Net weight: Approximately 250g			Parity: Even*, Odd, None
2. Ethernet interface			Stop bits: 1 bit* or 2 bits
Connector: RJ45		* : Factory settings	
Protocol: TCP/IP			

Weighing Instrument Cables


Required AD-8526 Cable by Weighing Instrument		
Weighing instrument	Interface option	Communications cable (Length 2 m)
GX, GF, GX-K, GF-K, GP, GR, HR, MC	None (D-Sub 25-pin, standard accessory)	AX-KO1710-200
EK-i, EW-i, FC-i, FC-S i, GH, HR-i	None (D-Sub 9-pin, standard accessory)	AX-KO2466-200
AD-4212C	None (D-Sub 9-pin, standard accessory)	AD4212C standard accessory
HV-G, HV-WP, HW-G, HW-WP	None (DIN 7-pin, standard accessory)	AX-KO1786-200
FG	OP-03 (DIN 7-pin)	AX-KO1786-200
FG-L, FG-M	OP-23 (DIN 8-pin)	AX-KO1786-200

Please refer to our website for the latest information.

Connecting the Weighing Instrument to a Computer

- Caution**
- Contact your network manager before connecting the AD-8526 to a network. The AD-8526 may generate a network error. A&D assumes no responsibility for any errors that may be caused.
 - Affix the AC adapter labels to the AC adapter and the connector end of the cable.
 - Set the IP address and subnet mask to the AD-8526 one at a time. The factory default IP address is 172.16.100.2. Do not duplicate an IP address.
 - The IP address cannot be returned to the factory setting. We recommend that you write the IP address on the supplied label.

- Preparation**
Please download each software from A&D website (<https://www.aandd.jp/products/software/software.html>) for below.
- Setting procedure "WinCT-Plus" instruction manual
 - IP address setting software "Device Installer"
 - Data acquisition software "WinCT-Plus"
- Connect the AC adapter to the AD-8526.
 - Connect the AD-8526 and the computer with a cross cable directly or connect them with a hub and straight cables.
 - Input the IP addresses and subnet masks of the AD-8526 and computer. For the setting procedure, refer to the "WinCT-Plus" instruction manual. **Communication error may occur if set incorrectly.**
 - Write the IP address on the IP address label and affix it on the AD-8526 (and the weighing instrument if necessary) where you can see the address easily.
 - Install the RsMulti data acquisition software to download the "WinCT-Plus" instruction manual.PDF form A&D website. (<https://www.aandd.jp/products/software/software.html>)
 - Connect the accessory RS-232C cable between the RS-232C interface of the weighing instrument and the AD-8526's port.

Data Communication Software

- The RSMulti data acquisition software can be launched on a Windows PC. Menu: [Start] → [Program] → [A&D WinCT-Plus] → [RsMulti]. Refer to the Manual PDF file in the A&D WinCT-Plus for details on the operation.
- This software can acquire data from multiple devices connected by LAN or RS-232C.
 - The software can control these devices with commands.
 - The software can acquire data transmitted from devices. Example: When the **PRINT** key is pressed on the balance, data is output and is acquired with the computer.
 - Stored data can be used with Excel. (Install Microsoft Excel before use.)

Example

A sample window of RsMulti

Two balances are connected.

	IP address	Subnet mask
Computer	172.16.100.1	255.255.0.0
GX-600	172.16.100.2	
GP-30K	172.16.100.3	